






KRISHZYME® Factor Xa Inhibitor Screening Assay Kit

REF : KBBA73

Ver1.0


RUO

Chromogenic Assay for Estimation of Factor Xa Inhibitor Activity in Human Blood Products

RUO	For Research Use Only	REF	Catalog Number
	Store At	LOT	Batch Code
	Manufactured By		Biological Risk
	Expiry Date		Consult Operating Instructions

For Research Use Only. Purchase does not include or carry the right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Krishgen Biosystems Private Limited is strictly prohibited.

REF KBBA73

 96 tests

Krishgen Biosystems Private Limited

For US / Europe: toll free +1(888)-970-0827 tel: +1(562)-568-5005
For Asia / India: tel: +91(22)-49198700
Email: sales1@krishgen.com

KRISHZYME® Factor Xa Inhibitor Screening Assay Kit

Introduction:

Factor X, also known as EC3.4.21.6 or thrombokinase, is a serine endopeptidase involved in blood coagulation. Factor X needs to be activated into Factor Xa to cleave prothrombin into active thrombin. Deficiency in inactive and active Factor X results in blood coagulation disorders. Warfarin, a common anticoagulant used in the treatment of thrombosis, inhibits the conversion of Factor X to Xa, resulting in lower risk of stroke and thrombosis.

Warfarin is used to treat pulmonary embolism and lower the risk of cardiovascular complications. The development of new Factor Xa inhibitors, also known as blood thinners, has been crucial in improving the quality of life of patients by bypassing the need for frequent blood tests.

Intended Use:

The KRISHZYME® Factor Xa Inhibitor Screening Assay Kit colorimetric assay designed to measure the activity of human Factor Xa for screening and profiling applications.

Principle:

The assay is based on a chromogenic method in which Factor Xa cleaves a specific synthetic substrate to release p-nitroaniline (pNA), producing a yellow colour measurable at 405 nm. To evaluate inhibitory effects, Factor Xa is preincubated with or without a test inhibitor prior to substrate addition. In the absence of an inhibitor, active Factor Xa cleaves the substrate efficiently, resulting in increased colour intensity, whereas in the presence of an inhibitor, enzymatic activity is reduced, leading to decreased pNA release and lower absorbance. Thus, the intensity of the colour is directly proportional to Factor Xa activity and inversely proportional to the inhibitory effect of the test compound.

Materials Provided:

1. Human Factor Xa (>1 ug) - 1 vial
2. Substrate – 50 ul
3. Dilution Buffer – 10 ml
4. 96 well plate - 1 Nos.
5. Instruction Manual

Materials to be provided by the End-User:

1. Microplate Reader / Spectrophotometer able to measure absorbance at 405nm.
2. Adjustable pipettes to measure volumes ranging from 10 ul to 2500 ul, duly calibrated.
3. 1.5 ml Eppendorf tubes.
4. Deionized (DI) water.
5. Software for data analysis, if required.
6. Timer/Stop watch.

Reagent Preparation, Storage, and Stability Information:

Human Factor Xa and the Substrate need to be stored in -80°C and store the Dilution Buffer at -20°C.

Note:

- Bring all reagents to room temperature.
- All reagents should be diluted immediately prior to use.

Preparation of test inhibitors:

- a. If the Test Inhibitor is water-soluble, prepare serial dilutions 10-fold more concentrated than the desired final concentrations in Dilution Buffer.
- b. If the Test inhibitor is soluble in DMSO, prepare the inhibitor in 100% DMSO at a concentration 100-fold higher than the highest desired concentration, then dilute the inhibitor 10-fold in Dilution Buffer to prepare the highest concentration of the 10-fold intermediate dilutions. The concentration of DMSO is now 10%. Use 10% DMSO in Dilution Buffer (vol/vol) as diluent for the serial dilution to keep the concentration of DMSO constant. For positive and negative controls, prepare 10% DMSO in Dilution Buffer (vol/vol) so that all wells contain the same amount of DMSO (Diluent Solution).

Note: The final concentration of DMSO should not exceed 1%.

Health Hazard Warnings:

1. The source materials have been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using FDA approved methods. However adequate care should be taken when handling these materials as a source for potentially infective agents.
2. The Buffer contains Sodium Azide as preservative. It may be harmful if ingested, inhaled or absorbed through the skin. Refer to product MSDS for details.
3. For Research Purpose and In-Vitro Laboratory Use Only.

Reagent Preparation:
1. Working Factor Xa Preparation:

- a. Thaw Factor Xa on ice. Briefly spin the tube to recover the full content
- b. Dilute Factor Xa to 0.125 ng/ul in Dilution Buffer (40 ul/well).

Note: Keep the diluted protein on ice until use. Discard any unused diluted protein after use.

2. Working Substrate Solution:

Dilute Substrate 100-fold in Dilution Buffer.

Assay Protocol:

Component	Negative Control	Positive Control	Test Inhibitor
Dilution Buffer	40 ul	-	-
Working Factor Xa	-	40 ul	40 ul
Test Inhibitor	-	-	10 ul
Dilution Buffer	10 ul	10 ul	-
30 minutes at Room Temperature			
Working Substrate	50 ul	50 ul	50 ul
Total	100 ul	100 ul	100 ul

1. After adding the Substrate, incubate the plate for 30 – 60 mins at room temperature OR perform Kinetic Analysis.
2. Take the reading at 405 nm.

Safety Precautions:

- **This kit is For Research Use only.** Follow the working instructions carefully.
- The expiration dates stated on the kit are to be observed. The same relates to the stability stated for reagents.
- Do not use or mix reagents from different lots.
- Do not use reagents from other manufacturers.
- Avoid time shift during pipetting of reagents.
- All reagents should be kept in the original shipping container.
- Source materials maybe derived from human body fluids or organs used in the preparation of this kit were tested and found negative for HBsAg and HIV as well as for HCV antibodies. However, no known test guarantees the absence of such viral agents. Therefore, handle all components and all patient samples as if potentially hazardous.
- Since the kit contains potentially hazardous materials, the following precautions should be observed.
 - Do not smoke, eat or drink while handling kit material.
 - Always use protective gloves.
 - Never pipette material by mouth.
 - Wipe up spills promptly, washing the affected surface thoroughly with a decontaminant.
- In any case GLP should be applied with all general and individual regulations to the use of this kit.



LIMITED WARRANTY

Krishgen Biosystems Private Limited does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the product; against defects in products or components not manufactured by Krishgen Biosystems Private Limited, or against damages resulting from such non-Krishgen Biosystems Private Limited made products or components. Krishgen Biosystems Private Limited passes on to customer the warranty it received (if any) from the maker thereof of such non-Krishgen made products or components. This warranty also does not apply to product to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Krishgen Biosystems Private Limited.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of Krishgen Biosystems Private Limited shall be to repair or replace the defective product in the manner and for the period provided above. Krishgen Biosystems Private Limited shall not have any other obligation with respect to the products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall Krishgen Biosystems Private Limited be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Krishgen Biosystems Private Limited with respect to the product. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

Krishgen Biosystems Private Limited. 2026

THANK YOU FOR USING A KRISHGEN PRODUCT!

KRISHGEN BIOSYSTEMS PRIVATE LIMITED®, DHARMAPLEX®, GENBULK®, GENLISA®, KRISHZYME®, KRISHGEN®, KRIBIOLISA®, KRISHPLEX®, TITANIUM®, QUALICHEK® are registered trademarks of KRISHGEN BIOSYSTEMS PRIVATE LIMITED.

©KRISHGEN BIOSYSTEMS PRIVATE LIMITED. ALL RIGHTS RESERVED.

KRISHGEN BIOSYSTEMS PRIVATE LIMITED | OUR REAGENTS | YOUR RESEARCH |