

KRISHZYME™ TEV Protease

Catalog Number: KBENZ71

Protein Description

The KRISHZYME™ TEV Protease (Tobacco Etch Virus nuclear-inclusion-a endopeptidase) is a highly sequence-specific cysteine protease from Tobacco Etch Virus (TEV). It is a member of the PA clan of chymotrypsin-like proteases. TEV Protease is a highly specific cysteine protease that recognizes the amino-acid sequence Glu-Asn-Leu-Tyr-Phe-Gln-(Gly/Ser) and cleaves between the Gln and Gly/Ser residues. Due to its high sequence specificity it is frequently used for the controlled cleavage of fusion proteins in vitro and in vivo.

Source:

Tobacco Etch Virus (TEV)

Expression Host:

E.coli

Purity:

>80% as determined by SDS-PAGE quantitative densitometry by Coomassie Blue Staining.

Endotoxin:

< 0.05 EU/1000 units as determined by the LAL method.

N terminal:

N- terminal Gst Tag

Molecular Mass:

The KRISHZYME™ TEV Protease has a calculated molecular mass of 50 kDa

Concentration:

5 U/ml

Unit Definition:

One Active Unit is defined as the reaction of 1X TEV Protease Buffer for 3 hours at 30°C and cleaves >35% of 3 ug substrate

Specific Activity:

17 U/ul

SDS-PAGE:

Fig.1.




Fig. 1. Purity analysis by SDS-PAGE Detection

Formulation:

KRISHZYME™ TEV Protease is supplied as a liquid containing 25 mM Tris-HCl, 50 mM NaCl, 1 mM EDTA, 50% glycerol. KRISHZYME™ TEV Protease is also supplied as a lyophilized enzyme. Please contact us by email at sales@krishgen.com for the same.

Reconstitution:

Being an enzyme, the concentration may differ from lot to lot. We always recommend referring the accompanying data sheet to view the exact concentration and the recommended dilution schemata.

Centrifuge the vial at 4°C before opening to recover the entire contents. Please contact us for any concerns or special requirements at +91-22-49198700 | Email: sales@krishgen.com

Storage:

Store it under sterile conditions at -20°C to -80°C upon receiving for at least 12 months. It is recommended to aliquot the enzyme into smaller quantities for optimal storage. Avoid repeated freeze-thaw cycles.

Application:

Removal of affinity purification tags such as maltose-binding protein (MBP) or poly-histidine from fusion proteins

Protocol:

Incubation times and enzyme concentrations must be optimized for a particular substrate. Typical reaction conditions are as follows:

- 1) Combine 15 ug of substrate and H₂O (if necessary) to make a 45 ul total reaction volume.
- 2) Add 5 ul of TEV Protease Reaction Buffer (10X) to make a 50 ul total reaction volume.
- 3) Add 1 ul of TEV Protease.
- 4) Incubate at 30°C for 1 hour or at 4°C overnight.

Notes:

1. In order to achieve the best digestion effect, please ensure that the recombinant protein is partially or completely purified protein.
2. If the fusion protein contains the denaturant, the denaturant should be removed to perform the digestion reaction.
3. To remove TEV Protease from the digested system, it is recommend to use a glutathione affinity column.
4. If the fusion protein sample contains >2 M urea, >0.5 M Guanidine hydrochloride, >50 mM imidazole, pH values below 6 or above 9, or cysteine protease inhibitors then it will be necessary to dialyze the fusion protein into TEV protease reaction buffer before TEV Protease cleavage.
5. TEV protease is inhibited by reaction buffers containing >40% Glycerol.
6. Inhibition occurs in the presence of ≥ 5 mM Zn^{2+} , ≥ 1 mM Cu^{2+} and ≥ 10 mM Co^{2+} .
7. Some substrates may require extended incubation periods (up to three days at either 4°C or 30°C) to achieve complete cleavage. The addition of more TEV Protease after 24 hours may also help achieve complete cleavage of some substrates.

References:

Improved solubility of TEV protease by directed evolution
S Van Den Berg, PÅ Löfdahl, T Hård... - Journal of biotechnology, 2006 - Elsevier

RB Kapust, DS Waugh - Protein expression and purification, 2000 - mcl1.ncifcrf.gov

Expression and Purification of Soluble His6-Tagged TEV Protease
JE Tropea, S Cherry, DS Waugh - High throughput protein expression and ..., 2009 - Springer

Cell-type-specific TEV protease cleavage reveals cohesin functions in Drosophila neurons
A Pauli, F Althoff, RA Oliveira, S Heidmann... - Developmental cell, 2008 - Elsevier

Enhancing the stability and solubility of TEV protease using in silico design
LD Cabrita, D Gilis, AL Robertson, Y Dehouck... - Protein ..., 2007 - Wiley Online Library

LIMITED WARRANTY

Krishgen Biosystems does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by Krishgen Biosystems, or against damages resulting from such non-Krishgen Biosystems made products or components. Krishgen Biosystems passes on to customer the warranty it received (if any) from the maker thereof of such non Krishgen made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Krishgen Biosystems.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of Krishgen Biosystems shall be to repair or replace the defective Products in the manner and for the period provided above. Krishgen Biosystems shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, and strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall Krishgen Biosystems be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Krishgen Biosystems with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

Krishgen Biosystems. 2025

THANK YOU FOR USING KRISHGEN PRODUCT!

KRISHGEN BIOSYSTEMS®, GENLISA®, DHARMAPLEX™, GENBULK™, GENLISA™, KRISHZYME®, KRISHGEN®, KRIBIOLISA®, KRISHPLEX®, TITANIUM®, QUALICHEK® are registered trademarks of KRISHGEN BIOSYSTEMS. ©KRISHGEN BIOSYSTEMS. ALL RIGHTS RESERVED.

KRISHGEN BIOSYSTEMS | OUR REAGENTS | YOUR RESEARCH |