



# Medizym<sup>®</sup> anti-CCP Ref

ELISA for the determination of IgG-antibodies against cyclic citrullinated peptides (CCP) in human serum and plasma



Early diagnosis of rheumatoid arthritis (RA)

Well proven CCP antigens of the 2. generation (peptides)

Calibrated with constant factor to WHO-standard W1066 (Rheumatoid factor)

Linearity in the range from 20 to 2000 U/ml

Sensitivity: 79 % (378 RA)

Specificity: 97 % (356 non-RA)

---

## KRISHGEN BioSystems