








KRISHGEN BioSystems



**Plant
Sciences**



**BT Cotton ELISA Kits
GMO ELISA Kits
Plant Antibodies
Tissue Culture Medias, Reagents
Plant Growth Mixes
Assay Components, Reagents
Conjugation Kits**



KRISHGEN BioSystems

Cat No	Product Particulars	Type of ELISA	Pack Size
Cry1Ab/Cry1Ac			
KBA001D	Cry 1Ab/Cry1Ac ELISA Deluxe	qualitative	1x96wells
KBA001D-10	Cry 1Ab/Cry1Ac ELISA Deluxe	qualitative	10x96wells
KBA001D-20	Cry 1Ab/Cry1Ac ELISA Deluxe	qualitative	20x96wells
KBA001S	Cry 1Ab/Cry1Ac ELISA Standard	qualitative	96wells
KBA001S-10	Cry 1Ab/Cry1Ac ELISA Standard	qualitative	10x96wells
KBA001S-20	Cry 1Ab/Cry1Ac ELISA Standard	qualitative	20x96wells
KBA001plate	Cry 1Ab/Cry1Ac Plate		1plate
KBA001conj	Cry 1Ab/Cry1Ac Conjugate		11ml
KBA001conj1	Cry 1Ab/Cry1Ac Conjugate		1Ltr
KBA002D	Cry 1Ab/Cry1Ac ELISA Deluxe	quantitative	96wells
KBA002D-10	Cry 1Ab/Cry1Ac ELISA Deluxe	quantitative	10x96wells
KBA002D-20	Cry 1Ab/Cry1Ac ELISA Deluxe	quantitative	20x96wells
KBA002S	Cry 1Ab/Cry1Ac ELISA Standard	quantitative	96wells
KBA001S-10	Cry 1Ab/Cry1Ac ELISA Standard	quantitative	10x96wells
KBA001S-20	Cry 1Ab/Cry1Ac ELISA Standard	quantitative	20x96wells
KBA002s	Cry 1Ab/Cry1Ac Standards Set		1x5x0.5ml
Cry 1Ab			
KBA003D	Cry 1Ab ELISA Deluxe	qualitative	1x96wells
KBA003D-10	Cry 1Ab ELISA Deluxe	qualitative	10x96wells
KBA003D-20	Cry 1Ab ELISA Deluxe	qualitative	20x96wells
KBA003S	Cry 1Ab ELISA Standard	qualitative	96wells
KBA003S-10	Cry 1Ab ELISA Standard	qualitative	10x96wells
KBA003S-20	Cry 1Ab/ ELISA Standard	qualitative	20x96wells
KBA003plate	Cry 1Ab Plate		1plate
KBA003conj	Cry 1Ab Conjugate		11ml
KBA003conj1	Cry 1Ab Conjugate		1Ltr
Cry1Ac			
KBA004D	Cry1Ac ELISA Deluxe	qualitative	1x96wells
KBA004D-10	Cry1Ac ELISA Deluxe	qualitative	10x96wells
KBA004D-20	Cry1Ac ELISA Deluxe	qualitative	20x96wells
KBA004S	Cry1Ac ELISA Standard	qualitative	96wells
KBA004S-10	Cry1Ac ELISA Standard	qualitative	10x96wells
KBA004S-20	Cry1Ac ELISA Standard	qualitative	20x96wells
KBA004plate	Cry1Ac Plate		1plate
KBA004conj	Cry1Ac Conjugate		11ml
KBA004conj1	Cry1Ac Conjugate		1Ltr
Cry2Ab			
KBA005D	Cry2Ab ELISA Deluxe	qualitative	1x96wells
KBA005D-10	Cry2Ab ELISA Deluxe	qualitative	10x96wells
KBA005D-20	Cry2Ab ELISA Deluxe	qualitative	20x96wells

KRISHGEN BioSystems

Cat No	Product Particulars	Type of ELISA	Pack Size
KBA005S	Cry2Ab ELISA Standard	qualitative	96wells
KBA005S-10	Cry2Ab ELISA Standard	qualitative	10x96wells
KBA005S-20	Cry2Ab ELISA Standard	qualitative	20x96wells
KBA005plate	Cry2Ab Plate		1plate
KBA005conj	Cry2Ab Conjugate		11ml
KBA005conj1	Cry2Ab Conjugate		1Ltr
Cry1c			
KBA006D	Cry1c ELISA Deluxe	qualitative	1x96wells
KBA006D-10	Cry1c ELISA Deluxe	qualitative	10x96wells
KBA006D-20	Cry1c ELISA Deluxe	qualitative	20x96wells
KBA006S	Cry1c ELISA Standard	qualitative	96wells
KBA006S-10	Cry1c ELISA Standard	qualitative	10x96wells
KBA006S-20	Cry1c ELISA Standard	qualitative	20x96wells
KBA006plate	Cry1c Plate		1plate
KBA006conj	Cry1c Conjugate		11ml
KBA006conj1	Cry1c Conjugate		1Ltr
Cry1F			
KBA007D	Cry1F ELISA Deluxe	qualitative	1x96wells
KBA007D-10	Cry1F ELISA Deluxe	qualitative	10x96wells
KBA007D-20	Cry1F ELISA Deluxe	qualitative	20x96wells
KBA007S	Cry1F ELISA Standard	qualitative	96wells
KBA007S-10	Cry1F ELISA Standard	qualitative	10x96wells
KBA007S-20	Cry1F ELISA Standard	qualitative	20x96wells
KBA007plate	Cry1F Plate		1plate
KBA007conj	Cry1F Conjugate		11ml
KBA007conj1	Cry1F Conjugate		1Ltr
Sundries			
KBAWB-01	Wash Buffer		120ml
KBAWB-02	Wash Buffer		1Ltr
KBAEXB-01	Extraction Buffer		500ml
KBAEXB-02	Extraction Buffer		1Ltr
KBASUB-01	Enzyme Substrate		11ml
KBASUB-02	Enzyme Substrate		1Ltr
KBASTOP	Stop Solution		11ml
KBASTOP	Stop Solution		1Ltr

KRISHGEN BioSystems

Note :

* *available
shortly*

** *ELISA Deluxe =
kit containing plate, conjugate, enzyme substrate, wash buffer & stop solution,
positive control or standards.*

*** *ELISA Standard =
kit containing plate, conjugate.*