



# Antiphospholipid Dot

## Profil Diagnostics for the Diagnose of Anti-Phospholipid Syndrom (APAS)

- Cardiolipin
- $\beta$ -Gylcoprotein I
- Phosphatidylserin
- Phosphatidylinositol

### Profil Diagnostics in Dot Format

- ✓ Simple and easy
- ✓ Especially convenient in case of low determinations
- ✓ No equipment needed
- ✓ No standard curve required
- ✓ Simultaneous screening IgG and IgM
- ✓ Extraordinary detection limit e.g. < 0.2 $\mu$ g/ml  $\beta$ GP I

Order Code: 5006

Interested? Please contact your local distributor:

**KRISHGEN BioSystems**