


ACTICLOT is a clot-based assay for determining resistance to activated Protein C (aPC) caused by the Factor V Leiden mutation (FV:Q506)

Attention Hematology Users – APCR ELISA kit now available!!!

Activated Protein C Resistance – ACTICLOT Assay

Protein C Anticoagulant Pathway



APCR ELISA

Samples: Citrate collected plasma
 Sample Preparation: Neat
 Size: 30 µL
 Assay Time: < 5 minutes
 Sensitivity: 100% for Factor V Leiden homozygous + heterozygous carriers

Precision: Intra-assay CV < 5% Inter- < 5%
 Assay CV
 No of Tests: 20
 Specificity: 100%
 Certifications: 510(k), CE
 Kit Components: PTA reagent (thrombin activator) RVV-V reagent
 RVV-V (+ APC) reagent
 Dilution plasma



Contact Us :
 135/37 Zaveri Bazar, Mumbai 400002.
 Tel: 91-22-66372990 Fax: 91-22-66372991
 Email: sales@krishgen.com <http://www.krishgen.com>