

## WASH BUFFER USED IN KRISHGEN ELISA KITS FAQs

ELISA wash buffer was developed as a high performance washing solution for use in all ELISA formats. The optimized formulation of pH stabilizers, salts and detergents increases the efficiency of washing stages during ELISA, leading to a reduction in background noise and an increased specific signal.

Store at +4°C.

DO NOT FREEZE

This product should be stored undiluted. Should this product contain a precipitate, warm the buffer in a 37°C water bath until dissolved. Do not boil.

Note: Cool to RT before using the same in the ELISA procedure.